## RAJIV GANDHI INSTITUTE OF PETROLEUM TECHNOLOGY, JAIS

DATE: 09th May 2021

## Corrigendum 2

Tender No: RGIPT/JAIS/CRF/E-OPN/LAB/2020-21/03 (Raman Spectrometer)

Please refer to the above e-tender no. RGIPT/JAIS/CRF/E-OPN/LAB/2020-21/04 on rgipt.euniwizarde.com portal & institution website. The dates of the bidding process for GC-FID-TCD-DHA is revised as follows:

- The revised technical specifications of the above is uploaded in the website
- The revised date for the last date of submission of e-tender is 19th June 2021, 4 PM.
- The date of opening of technical bid is 19th June 2021, 5 PM

Other terms and conditions remain unchanged.

Note: Please see www.rgipt.ac.in for more information in the bid.

U. Ojha

Chairman, CRF Purchase Committee

Technical Specifications	Bidders Specificatio n	Complianc
Raman System		
250 mm or more focal length spectrophotometer		
Spectral Range- 50- 4000 cm 1 or more, Wave length		
range 200 - 2000 nm		
Spectral Resolution - 0.5 cm or better FWHM Appropriate inbuilt calibration source at least two, for		
		1 - 10 1 1 1 1 1 1
full spectral range. Spectrophotometer with high throughput (Kindly		
mention the % transmittance value)		
Gratings: 2400 gr/mm, 600 gr/mm ruled gratings		
Photoluminescence should be able to perform. PL range		
500-1050 nm		
Objectives: 10x, 100x, 50x Raman filter set for 532 nm, ~50 cm   edge filter, beam		
Kaman niter set for 552 mil, ~50 cm - edge meet, boats		
splitter and band pass filter  ND filter or alternate technique to vary		
laser power 1 to 100%.	577 6 31 A. S	
Camera for sample viewing: ~3 Mega Pixel or better		
Laser		
Laser 532 nm, Diode Laser, 50 mW power or Higher at		
laser output		1
Laser power stability ≤ 1% or Less		E
Beam Diameter ≤ 1.2 mm		
Laser Line Width < 0.00019 nm or better		
100: 1 Polarization stability	E STATE OF THE STA	
Power control must be software controlled preferably		and the second second
Fiber coupled or Direct coupled to microscope		
Detector		dans gas-n
CMOS or Back illuminated CCD Detector		
1 stage TE cooled to -60 °C or better		
QE > 90% @530 nm		
Pixel Size: 26 µm x 26 µm		
Dark Noise 0.001 e-/pixel/s		-
USB Interface		
Software		
System- Window 10 based or Higher- License window 10 required		1
Data acquisition software for instrument control, image		
capture for viewing and saving of images.		
Data processing software for baseline correction, auto		
fluorescence correction, sensitivity correction, curve		

1915)21 1915)21

Penas Jasta

Harzana Al 19-8-102)

smoothing, etc., with other compatible application	
Accessories	
One Computer with the latest configuration- Dual core Processor or better, 8GB RAM, 1000 GB HDD or better:  24" or Higher LCD monitor, suitable UPS system	
and Wibration Table: Closed champer	
(Microscope/detector in closed system)	and a second sec
Warranty	
3 years from the date of successful installation	The state of the s
Installation and Training On-site Installation and Training by expert application	
Engineer  Equivalent system should have been supplied to well- known research institution including IITs, IISc and CSIR	
Note: The system should have facility to upgrade in	
future as per following requirements; 325 nm He-Cd laser, 785 nm Diode laser; motorized scanning stage with integrated measuring system; Microscope heating cooling stage; or any other (provide complete user list of such items).	

Ampan Jal

M3. 2 (195/19

Harry A Zariza

Shila .. 21

Um